


# TICK CARD

## INSTRUCTIONS FOR USE


### IMPORTANT!

- Remove ticks as soon as possible.
- Do not squeeze, twist or burn a tick.
- Do not use petroleum jelly, oil or any other chemical products to remove ticks.
- Stressing a tick greatly increases the risk of it injecting harmful agents into the bloodstream of its host.




- 1.** Slip the notch of the Tick Remover Card under the tick.

NOTE: For smaller ticks or ticks that are not yet fully engorged, use the smaller notch in the Tick Remover Card.



- 2.** Slowly push the card forwards and upwards with gentle and even pressure.

NOTE: Do not apply excessive force as this may stress the tick, causing the head to remain in the skin.



- 3.** Allow the tick to let go and lift away from the skin with ease.

NOTE: Disinfect the bite area and the Tick Remover Card after removal using hygienic cleansing wipes.

### TICK-BORNE DISEASE

A small red itchy patch is a normal body-reaction to a tick bite but symptoms can be much more severe. Serious tick-borne diseases include Lyme disease (Lyme borreliosis) which is increasingly common throughout Europe and America.

### SYMPTOMS

Symptom patterns vary from person to person. Early indications in the days after the bite can include flu-like symptoms, extreme tiredness, headache, muscle pain, tender glands and sensitivity to temperature, sound and light. In some cases you may see a roughly circular red rash, without itching or pain, that slowly spreads outwards. It may be clear in the middle and is often called a bull's eye rash. Seek medical advice promptly if you are concerned or experience any of the above symptoms after a bite.